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*For English see other side

Instructions STD test men (urine)

With this test kit you can test yourself for chlamydia (urine). Check before you test if the kit is complete. Also check the expiration date on the test tube and the packaging of the cotton swabs.

The chlamydia test for men (urine) contains:

- A test tube filled with transport liquid
- A pipet
- An application form
- A white (plastic) shipping envelope
- A Rigid Safety Bag
- A tube holder

Instructions:

1. Get the application form and write at 'afnamedatum' which date it is today (the date you're carrying out the test).
2. Take a **clean** cup and collect your first drops/stream of urine in this cup (max. 30 ml).
3. Remove the cap of the test tube by rotating it, and place the tube in the tube holder. Make sure that the liquid stays inside the tube.
4. Use the pipet to collect a part of the urine, and proceed by putting it in the test tube until the window is filled halfway through (so fill until it is between the two black lines).
5. Close the test tube by firmly rotating the cap onto the tube.

6. Place the test tube according to the instructions in the Rigid Safety Bag. Remove the white strip on the Rigid Safety Bag to seal the bag.
7. Put the Rigid Safety Bag in the white (plastic) shipping envelope.
8. **Please note: also include the application form in the shipping envelope!**
9. Send the shipping envelope (including the Rigid Safety Bag with the test tube, and the application form) to the lab by placing it in an orange PostNL postbox before 5 pm on a weekday. You don't have to put stamps on it.
 - Is the outside temperature below 0 or above 30 degrees celsius? In that case, bring the package to a PostNL postal office.
10. After 3-4 business days go to the results page on [SoaStudentArts.nl](https://www.soastudentarts.nl). Fill in your client number and you will receive your results anonymously.

Contact

E-mail: info@soastudentarts.nl

Tel: +31 (0) 50 211 0 911

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